

Central Vaccination Clinics: TrackVax V 23 New Functionality: Syringe Mode

There is now an option to supply more than one vial at a time in syringe mode.

Process Steps

Follow Instructions on-screen:

1. After scanning the box to supply trays and syringes from, there will be a new option, to select the number of vials required. (Similar to the functionality in Vial Mode.)

The screenshot shows the TrackVax Vial Labelling Station interface. At the top, it displays 'Print Tray & Syringe Labels . . .', '27-Oct-2021 10:14', and a 'Clear' button. Below this, there are fields for 'Labelled By:', 'Checked By:', and 'Current Batch:'. A 'Back' button is also present.

The main interface is divided into several sections:

- Product Information:** SSSC (00): 153915394701781143, GTIN (01/02): 04260703260002, COVID-19 mRNA Vaccine BNT162b2 30micrograms/0.3ml do, Use By (17/7003): 31-Oct-2021 15:16, Batch (10): FG9019, Serial No (21):, Reference (251): 4004194/6, Vial No. (250):, Discard By (7003): 27-Oct-2021 16:13.
- Label:** Comirnaty Pfizer-BioNTech, *** DILUTED COVID-19 VACCINE VIAL ***, *** TEST *** PF, TRAY NO., No. of Syringes: 0, (0104260703260002), (70030110071610), (17000010), Vial No:, Batch: FG9019, Discard By: 27-OCT-2021, *** TEST *** 16:13.
- Saline GTIN:** 04030539068079, **Description:** 0.9% w/v Sodium Chloride Injection, **Batch Number:** 00077015, **Expiry Date:** 31-Aug-2022.
- Issue To Bay:** Mode, **No. of Vials (Max 12):** 2, **Split Trays**, **Print 2 Tray Label & 14 Syringe Labels**.

Instructions: Select the No. of Vials required and then Scan the 1st Staff Member. [Click here to scan source label . . .](#)

TrackVax Main Menu | 2300 | C8F7503C895A | Logout in 30 Mins | 27-Oct-2021 10:14 | Ver: 23

This close-up shows the 'No. of Vials (Max 12):' section of the interface. It features a numeric input field with the value '2', flanked by left and right arrow buttons. To the right of the input field are two buttons: 'Split Trays' and 'Print 2 Tray Label & 14 Syringe Labels'.

Central Vaccination Clinics: TrackVax V 23 New Functionality: Syringe Mode

There is now an option to supply more than one vial at a time in syringe mode.

Process Steps

Follow Instructions on-screen:

- Once the number of Vials is selected (maximum will be the number available in the selected box or 12), the option to split trays remains.

The screenshot displays the TrackVax software interface. On the left, a 'Split Vaccine Trays' dialog box is open, showing a grid of 16 buttons numbered 1 to 16. Below the grid, it indicates '3 Trays' and '2 + 2 + 3'. An 'Ok' button is visible. On the right, a vaccine label is shown for 'Comirnaty Pfizer-BioNTech'. The label includes a QR code, 'TRAY NO.', 'No. of Syringes: 0', 'Vial No.', 'Batch: 20377015', and 'Discard By: 27-OCT-2021 21:07'. Below the label, there are fields for 'Batch Number: 20377015' and 'Expiry Date: 31-Aug-2023'. At the bottom, there are buttons for 'Split Trays' and 'Print 6 Tray Labels & 14 Syringe Labels'. The interface also shows a 'Current Batch' field and a 'Clear' button at the top right.

- This split will apply to all vials being supplied in this action; e.g. If issuing 2 vials with 3 trays per vial of $2+2+3=6$ tray labels will be printed with 14 syringe labels.
- There will be a gap of a blank label between each set of tray and syringe labels which come from the printer, and the vial number will remain visible on the label- apply the labels to the trays and syringes, in order of the corresponding vials.

For any questions, please contact your local PFH engineer who can raise a ticket using Service Cloud.