

There is new functionality in TrackVax to cater for the 0.5ml “Standard Dose” and 0.25ml “Booster Dose”


### Process Steps

***Before issuing Moderna in a session please confirm that you are set to the correct Mode of operation for TrackVax.***

1. Go to Stock Management from the main menu



2. On the Moderna Panel (Moderna must be in stock) a new button will appear to allow a Moderna vials to be issued in either “Booster Mode” or “Standard Mode”.
3. To change between modes- click on this button and scan two active staff badges.



There is new functionality in TrackVax to cater for the 0.25ml “Booster Dose” and 0.5ml “Standard Dose”

## Process Steps

**Before issuing Moderna in a session please confirm that you are set to the correct Mode of operation for TrackVax.**

- Proceed to supply, issue and return vials as normal in vial or syringe mode.
- When a Moderna box is scanned to supply labelled vials or syringes and trays it will clearly state that you are operating in Booster Mode when in that mode - see below:

The screenshot shows the 'Print Vial Labels' interface. At the top, it displays 'Print Vial Labels . . .', '24-Nov-2021 18:48', and a 'Clear' button. Below this are fields for 'Labelled By:', 'Checked By:', and 'Current Batch:', each with a redacted area. A '< Back' button is on the left. The main area is divided into two columns. The left column contains fields for 'SSCC (00):', 'GTIN (01/02): 30380777700688', 'Moderna COVID-19 VACCINE \*\* BOOSTER \*\*', 'Use By (17/7003): 12-Dec-2021 23:59', 'Batch (10): 46TEST', 'Serial No (21): 5588999RREEE', 'Reference (251):', 'Vial No. (250):', and 'Discard By (7003): 24-Nov-2021 23:59'. The right column shows a 'Label' preview for 'Moderna \*\* BOOSTER \*\*' (MB) with a QR code and various details like 'Labeled: 24-NOV-2021 18:48', 'Batch Number: 46TEST', and 'Discard By: 24-NOV-2021 23:59'. Below the label preview is a purple banner with '\*\*\* BOOSTER MODE \*\*\*'. At the bottom, there's a 'Mode' selector set to 'Booster Mode', 'No. of Vials (Max 2): 1', and a 'Print 1 Vial Label' button. A red banner at the very bottom contains instructions: 'Instructions: Select the No. of Vials required and then Scan the 1st Staff Member Badge . . . Click here to scan source label . . . Select Box'.

- Depending on the mode you are in the label will contain:
  - “**Standard Mode**”: Moderna Standard (MO) and 0.5ml per dose
  - “**Booster Mode**” : Moderna Booster (MB) and 0.25ml per dose

There is new functionality in TrackVax to cater for the 0.25ml “Booster Dose” and 0.5ml “Standard Dose”

### Process Steps

7. When Moderna is supplied in “Booster Mode” it will appear on the Dashboard in a new ‘light purple’ panel as Moderna Booster.

	Comimaty PF		Moderna BOOSTER	
	<input type="button" value="Syringes"/>	<input type="button" value="Vials"/>	<input type="button" value="Syringes"/>	<input type="button" value="Vials"/>
Supplied:	11	1	76	1
Returned:				
In Use:	11	1	76	1
Expires In:				
Expires In:				
Out of Date:				
Avg. Vial Yield:				
Actual Yield:				
Discarded:				
	* 22-Nov-2021 (-2 Days)*			

 Refresh (5 Mins)

8. Moderna Stock will only appear as in stock under the normal Moderna Panel- there will be **no stock in the Moderna Booster panel under stock management**
9. When returning Moderna Booster Doses the number of doses will default to 24. (this is subject to CVC Medication Management Guidance which should always be referred to.)
10. Moderna Booster doses will be identifiable in usage reports and in End of Session reporting as Moderna Booster.

*For any questions, please contact your local PFH engineer who can raise a ticket using Service Cloud.*