

Name of Test	Referral Lab	Sample Type	Special Requirements
ACTH	St James	EDTA x 2	< 4hr Frozen
C peptide	St James	Serum	< 1hr Frozen
Aldosterone	St James	EDTA x 2	<1hr Frozen
GAD antibody	St James	Serum	< 1hr Frozen
Gastrin	St James	Serum	< 4hr Frozen
IgF-1	St James	Serum	< 1hr Frozen
Insulin	St James	Serum	< 1hr Frozen
Intrinsic Factor	St James	Serum	< 4hr Frozen
PTH	St James	Serum	< 4hr Frozen
Renin	St James	EDTA x 2	<1hr Frozen (Morning preferred)
Alkaline Phosphatase– Bone Specific	Biomnis	Serum	<1hr Frozen
Amikacin	Biomnis	Serum	<4hr Frozen
Amino acids	Biomnis	1ml Lithium Heparin	<1hr Frozen
Anafranil (Clomipramine)	Biomnis	Serum	<1hr Frozen
Anti-Acetylcholine Receptor antibodies	Biomnis	Serum	<1hr Frozen
Antidiuretic Hormone (Vasopressin)	Biomnis	EDTA	<1hr Frozen
Chromogranin A	Biomnis	Serum	<1hr Frozen
CTX/ Crosslaps	Biomnis	Serum	<4 hr Frozen
Ethylene Glycol	Biomnis	Serum	<4 hr Frozen
Free Testosterone	Biomnis	Serum	<4 hr Frozen
Homocysteine Level	Biomnis	Serum	<1hr Frozen
Keppra (Levetiracetam)	Biomnis	Serum	<1hr Frozen
Lamictal (Lamotrigine)	Biomnis	Serum	<1hr Frozen
Lipoprotein A	Biomnis	Serum	<4 hr Frozen
Metanephrines	Biomnis	2 X 2ml Lithium Heparin	<1hr Frozen
Neurontin (Gabapentine)	Biomnis	Serum	<4 hr Frozen
Procollagen III	Biomnis	Serum	<1hr Frozen
Vasoactive Peptide	Biomnis	EDTA + aprotinine tube	<4 hr Frozen
Vitamin A	Biomnis	2ml Heparin plasma	<1 hr Frozen
Vitamin B1 (Thiamin)	Biomnis	EDTA protected from light	<4 hr Frozen
Vitamin B2 (Riboflavin)	Biomnis	EDTA protected from light	<4 hr Frozen
Vitamin B6	Biomnis	EDTA protected from light	<4 hr Frozen
Vitamin C	Biomnis	3ml Lithium Heparin protected from light	< 1hr Frozen
Bone Bank Serology	NVRL	Serum	Frozen- if received after 1pm